

Supplies for Specimen Collection and Transport*

BETTER DIAGNOSTICS.
BETTER PATIENT CARE.

APMG offers a variety of products to collect and transport the most common test specimens.



Thin Prep Pap Test

- Tests:
- Thin Prep Pap
 - HPV
 - Chlamydia/Gonorrhea

Stable at room temperature for up to 30 days after collection. Expiration date on vial.
DO NOT REFRIGERATE.



Conventional Pap Slide Kit

- Test: ● Conventional Pap
- Stable at room temperature once slide is fixed. No expiration date.
DO NOT REFRIGERATE.




Digene DNA Collection

- Tests: ● HPV (high and low)
● Chlamydia/Gonorrhea
- Stable at room temperature for up to 30 days after collection. Expiration date on package.
DO NOT REFRIGERATE.



Vaginosis Panel (BD Affirm VP III)

- Tests: ● Candida species
● Gardnerella vaginalis
● Trichomonas vaginalis
- Stable at room temperature for up to 3 days after collection. Expiration date on package.
DO NOT REFRIGERATE.



APMG Swab Collection Kit

- Tests: ● Cystic Fibrosis
● Group B Streptococcus
● Herpes Simplex 1&2
- Stable at room temperature for up to 30 days after collection. Expiration date on vial.
DO NOT REFRIGERATE.

Supply Collection Containers




Thin Prep Pap Test Vial with Broom



Thin Prep Pap Test Vial with Brush and Spatula



10% Formalin Pre-filled Tissue Biopsy Containers



APMG

Associated Pathology
Medical Group, Inc.

800.848.2764

apmglab.com

*More products are available. Call APMG for information.