

# Supplies for Specimen Collection and Transport\*

BETTER DIAGNOSTICS.  
BETTER PATIENT CARE.

APMG offers a variety of products to collect and transport the most common test specimens.



## Thin Prep Pap Test

- Tests:
- Thin Prep Pap
  - HPV
  - Chlamydia/Gonorrhea

Stable at room temperature for up to 30 days after collection. Expiration date on vial.  
DO NOT REFRIGERATE.



## Conventional Pap Slide Kit

- Test:
- Conventional Pap
- Stable at room temperature once slide is fixed. No expiration date.  
DO NOT REFRIGERATE.




## Digene DNA Collection

- Tests:
- HPV (high and low)
  - Chlamydia/Gonorrhea
- Stable at room temperature for up to 30 days after collection. Expiration date on package.  
DO NOT REFRIGERATE.



## Vaginosis Panel (BD Affirm VP III)

- Tests:
- Candida species
  - Gardnerella vaginalis
  - Trichomonas vaginalis
- Stable at room temperature for up to 3 days after collection. Expiration date on package.  
DO NOT REFRIGERATE.



## APMG Swab Collection Kit

- Tests:
- Cystic Fibrosis
  - Group B Streptococcus
- Stable at room temperature for up to 30 days after collection. Expiration date on vial.  
DO NOT REFRIGERATE.

## Supply Collection Containers




Thin Prep Pap Test Vial with Broom



Thin Prep Pap Test Vial with Brush and Spatula



10% Formalin Pre-filled Tissue Biopsy Containers



**APMG**

Associated Pathology  
Medical Group, Inc.

800.848.2764

apmglab.com

\*More products are available. Call APMG for information.